

Tilt-Pro™



OIS

Oncology Imaging Systems

Tilt-Pro™ TILTING BASE

The Tilt-Pro Tilting Base allows precise angling of the head in both prone and supine treatments to avoid critical structures. The Tilt-Pro adapts easily to most treatment and CT couch tops. It is designed for one-step molding, positioning, and indexing of Aquaplast RT™ and Fibreplast™ U-Frame masks.

Indexable channel arms allow adjustments from 0 to 30 degrees, in 5 degree increments. The Tilt-Pro Tilting Base uses traditional or Slimline U-Frames, Silverman Head Supports, and MOLDCARE® Head Cushions.



MR Tilt-Pro™ TILTING BASE

The Tilt-Pro is available in an MR compatible version to enable MR imaging in the same position as treatment.

PRONE CSA KIT


The Prone CSA Kit greatly simplifies the molding of prone masks. The kit includes everything needed to use your Tilt-Pro Base for Prone Cranial Spinal Axis treatment, including a molding fixture (for forming the prone mask), prone mask foam (to cushion the prone mask), and Double-High Nylon Swivel Locks (for securing two masks at once).



ProBoard™ FLAT HEADBOARD

The Tilt-Pro Base's carbon fibre top plate (the ProBoard) is available as a standalone device for standard flat treatments.





OIS is the distributor for Qfix in the UK & Ireland.

TiltPro, ProBoard, Aquaplast RT and Fibreplast are trademarks of Qfix Systems. MOLDCARE is a trademark of ALCARE CO, LTD.

DOCUMENT NO. BR52 / REV 2008 | For the latest information, please visit our website.



+44 (0) 1825 840633
sales@oncologyimaging.com
www.oncologyimaging.com