

**OIS launch new quality control software for radiotherapy and medical imaging  
QUALIMAGIQ® from QualifFormD SARL**

Developed in association with the Image Video and Communication Team at the *Research Institute of Communication and Cybernetics* in Nantes (IRCCyN), QUALIMAGIQ® is a modular software solution which allows you to organise and analyse the quality control procedures (QC) of Radiotherapy Treatment devices, Computed Tomography scanners, Gamma-Cameras, Digital Radiology, and MRI devices.



QUALIMAGIQ® is easy, comprehensive and precise. The software is built around a common platform where interfaces, methods of analysis and operating concepts are identical, regardless of the modality. It has 2 modes for the recovery, decoding and storage of DICOM images. The *Automatic Mode* takes roughly 10 seconds to perform; you click to choose the desired images, click again to analyse, then one last time to validate the analysis and generate 2 .pdf reports (in comprehensive or simplified mode), create a statistics point and transfer the images from the work file to an archive file. The step-by-step *Manual Mode* allows you to customise measurement masks and calculation methods, it's also a great introduction to QA procedures and makes a wonderful training tool.

QUALIMAGIQ® recognizes 3 user levels: consultation only, consultation with validation/rejection options and full user capacity, that allows you to regulate protocols, analysis methods, target values and tolerances. After sales service offers personalised off-site regulation via a customer parameters file and ftp account, so you can upload a problem image that is regulated and sent back for reintegration. QUALIMAGIQ® software comes with full and customer specific training.

If you'd like to see QUALIMAGIQ® in action at their UK product launch, come to the Oncology Imaging Systems Ltd. booth at the 2010 IPEM Biennial Radiotherapy Meeting, in the Julian Hodge Building, Cardiff University on the 6th-7th July. You'll be able to discuss the software with the engineers who designed it and get a full demo of all functions.

**For more information on QUALIMAGIQ® and all our patient positioning and immobilisation devices, please contact Oncology Imaging Systems Ltd., tel. 01825 840 633 fax. 01825 841 496 [sales@oncologyimaging.com](mailto:sales@oncologyimaging.com) [www.oncologyimaging.com](http://www.oncologyimaging.com)**