

Mould Room Couch



OIS

Oncology Imaging Systems

Qfix MOULD ROOM COUCH

The Mould Room couch combines an electronic, height-adjustable base on lockable castors with a fully indexed carbon fibre top. This easy-to-use and mobile couch allows patient set-up techniques to be mirrored in a Mould Room environment, without restrictive treatment or simulation room access. The rigid carbon fibre top is resistant to damage from water, plaster, other moulding products or bodily fluids and the whole unit can be lowered to a minimum height of 43cm, perfect for getting patients on and off the couch.

The Mould Room Couch base comes with locating bolts so it can be placed in the mould room at a fixed position to work in conjunction with lasers, which gives a comparable set up procedure in all areas (mould room, CT, linac); or it can be wheeled around to suit your needs.



New: Battery Powered Option

New battery-powered option includes two batteries: one to be charged while the other is in use. Ideal for room situations where power cabling is complex or may create a trip hazard.





OIS is the distributor for Qfix in the UK & Ireland.

DOCUMENT NO. BR50 / REV 2012 | For the latest information, please visit our website.



+44 (0) 1825 840633
sales@oncologyimaging.com
www.oncologyimaging.com